



	Model type
	LT
H (mm)	631
B (mm)	1100
T (mm)	594
T1 (mm)	1119
Doors	1
Tare weight (kg)	195
Inside cabinet volume (m ³)	0,11

* Standard base height (see page 4)

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: COOL UTS ergo LT Type 90, safety cabinet with wing doors		
	drawing no.: 000107636	page: 1 of 5	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	date: 2020-09-03		




	Model type
	LT
H1 (mm)	266
BS (mm)	384
TS (mm)	338
M (kg)	20
Max. additional load (kg)	20
L (l)	4
V (l)	5,3
LB (mm)	388

H1: Max. setting height
 BS: Width storage shelf
 TS: Depth storage shelf
 M: Load bearing capacity (uniformly distributed)
 L: Max. volume of the individual container
 V: Max. collection volume
 LB: Clear width
 TE: Technical unit:
 Climate class SN
 IP code IP21
 Refrigerant R 134a
 Temperature range: +2°C bis +10°C
 Supply voltage 230V / supply frequency 50Hz
 Power supply cable length approx. 2m, plug type E+F, CEE 7/7

Inner dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: COOL UTS ergo LT Type 90, safety cabinet with wing doors		
	drawing no.: 000107636	page: 2 of 5	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2020-09-03	




NW: 75
 AB: Exhaust air connection
 ZU: Supply air




Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: COOL UTS ergo LT Type 90, safety cabinet with wing doors		
	drawing no.: 000107636	page: 3 of 5	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein			

Classical base



ToeKick-base




Variable base versions		
	Classical base	ToeKick-base
HS (mm)	20 bis 135	60 bis 135


Base versions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: COOL UTS ergo LT Type 90, safety cabinet with wing doors		
	drawing no.: 000107636	page: 4 of 5	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 00	date: 2020-09-03	

Ex-zone 2
in cooler compartment



Ex-zone 2
in cooler compartment



Ex-zones

This drawing is our property!
Furnish this information to other
persons and reproduce or copy is
not permitted without our
written permission.

title: COOL UTS ergo LT
Type 90, safety cabinet with wing doors

drawing no.: 000107636 page: 5 of 5

revision: 00 date: 2020-09-03

Düperthal Sicherheitstechnik
Frankenstrasse 3
D-63791 Karlstein

 **DÜPERTHAL**

