







	Model type		
	M	L	XL
H (mm)	2114	2114	2114
B (mm)	594	894	1194
T (mm)	612	612	612
T1 (mm)	1112	1018	1168
Doors	1	2	2
Tare weight (kg)	69	91	113
Inside cabinet volume (m³)	0,7	1,0	1,3

Outer dimensions

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: Safety cabinet CHEM standard und CHEM tox		
	drawing no: D00104471	page: 1 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2020-02-25	




	Model type		
	M	L	XL
BS (mm)	584	884	1184
TS (mm)	560	560	560
M (kg)	50	50	50
Max. additional load (kg)	200	200	200
L (l)	20	30	40
V (l)	26	34	46
LB (mm)	494	765	1065
Number of storage shelves	3	3	3

BS: Width storage shelf
 TS: Depth storage shelf
 M: Load bearing capacity per storage shelf (uniformly distributed)
 L: Max. volume of the individual container
 V: Max. collection volume
 LB: Clear width

Max. container height (mm)		
		Model type
		M, L, XL
Storage shelves	H1	600
	H2	600
	H3	552
Storage shelves adjustable in 25mm pitch		

Inner dimensions


This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: Safety cabinet CHEM standard und CHEM tox		
	drawing no.: 000104471	page: 2 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2020-02-25	



NW: 110, with adapter 75
AB: Exhaust air connection

Ventilation connection

This drawing is our property! Furnish this information to other persons and reproduce or copy is not permitted without our written permission.	title: Safety cabinet CHEM standard und CHEM tox		
	drawing no: 000104471	page: 3 of 4	
Düperthal Sicherheitstechnik Frankenstrasse 3 D-63791 Karlstein	revision: 01	date: 2020-02-25	



● Loaded footprint

Geometry of the adjustable feet

This drawing is our property!
Furnish this information to other
persons and reproduce or copy is
not permitted without our
written permission.

Düperthal Sicherheitstechnik
Frankenstrasse 3
D-63791 Karlstein


title: Safety cabinet
CHEM standard und CHEM tox

drawing no.: 000104471

page: 4 of 4

revision: 01

date: 2020-02-25



 **DÜPERTHAL**