




	Safety cabinet model type SUPREME		
	L	XL	XXL
Article number storage shelf	79-80960-122	79-81260-122	79-81460-122

This drawing is our property!
Furnish this information to other
persons and reproduce or copy is
not permitted without our
written permission.

Düperthal Sicherheitstechnik
Frankenstrasse 3
D-63791 Karlstein


title: Installation instruction storage shelf
SUPREME L, XL, XXL

drawing no.: D00110402





page: 1 of 2

revision: 02

date: 2021-02-17



 **DÜPERTHAL**



Safety cabinet model type SUPREME			
	L	XL	XXL
A1 (mm)	638,5	938,5	1138,5

Parts list		
Pos.	Title	Quantity
010	Storage shelf	1
020	Bracket	1 pair
030	U-rail	1 pair
040	Chipboard screw countersunk head 4x40	12
050	Cylinder screw M4x10	4
060	Hexagon nut M4	4
070	Knurled nut	2
080	Chipboard screw half rounded head 4x25	2

This drawing is our property!
Furnish this information to other
persons and reproduce or copy is
not permitted without our
written permission.

Düperthal Sicherheitstechnik
Frankenstrasse 3
D-63791 Karlstein


title: Installation instruction storage shelf
SUPREME L, XL, XXL

drawing no.: D00110402

page: 2 of 2

revision: 02

date: 2021-02-17



DÜPERTHAL