



Product Service

CERTIFICATE

No. Z1A 012906 0510 Rev. 00

Holder of Certificate: **DÜPERTHAL SICHERHEITSTECHNIK
GMBH & Co. KG**



DÜPERTHAL

Frankenstraße 3
63791 Karlstein
GERMANY

Certification Mark:



Product: **Safety cabinets**

Tested according to: DIN EN 14470-1:2004
EK5/AK4 09-10:2009
DIN EN 16121:2017
DIN EN 16122:2012
AfPS GS 2019:01 PAK
EK5/13-11

The product meets the safety and health requirements of the German Product Safety Act section 20 to 22 ProdSG. The certification marks shown above can be affixed on the product. It is not permitted to alter the certification marks in any way. In addition, the certificate holder must not transfer the certificate to third parties. This certificate is valid until the listed date, unless it is cancelled earlier. All applicable requirements of the testing and certification regulations of TÜV SÜD Group have to be complied. For details see: www.tuvsud.com/ps-cert

Test report no.: 713188595

Valid until: 2025-08-12

Date, 2020-09-28

(Frank Feihle)

CERTIFICATE

No. Z1A 012906 0510 Rev. 00

Factories: 062099

Model(s): UTS ergo, BENCH

Parameters:

UTS ergo	
model type	approx. size without socket in mm (width x depth x height)
S and ST	601 x 594 x 596
M and MT	889 x 594 x 596
L, LT and LD	1100 x 594 x 596
XL and XLT	1400 x 594 x 596
S-5 and ST-5	601 x 494 x 596
M-5 and MT-5	889 x 494 x 596
L-5, LT-5 and LD-5	1100 x 494 x 596
XL-5 and XLT-5	1400 x 494 x 596

BENCH	
model type	approx. size without socket in mm (width x depth x height)
S and ST	601 x 594 x 770

Class of fire resistance: FWF 90

A detailed description of the specifications can be found in the test report.